



biolight 7000sc

Professional biolight 7000sc USER MANUAL

Professional Light Therapy system

HIGH OUTPUT, COLOR SELECTION
PULSE SELECTION, INTENSITY SELECTION

Revitalize Your Radiance with biolight 7000sc:
Selections of color setting

1. Red
2. Blue
3. Yellow
4. Green
5. Red & Blue
6. Red & Yellow
7. Red, Blue, Yellow & Green



RICHWAY & FUJI BIO INC.

1314 SOUTH KING STREET SUITE 520, HONOLULU, HI 96814 USA

TEL: +1-808-589-2800 TOLL-FREE: +1-855-338-6410 FAX: +1-808-597-1651

www.richwayandfujibio.com



 **biolight** 7000sc

Experience with Advanced LED Technology

Discover the biolight 7000sc and unlock the potential of advanced LED technology for your skincare journey.

Unlike ordinary LED light therapy devices that offer weak output and limited functionality, biolight 7000sc has been developed as a high-quality skincare device. It delivers comparable performance to professional-grade equipment used in specialized clinics.

Here's what sets biolight 7000sc apart

1] High Output for Enhanced Skincare:

Experience the power of strong output, significantly boosting the effectiveness of your skincare routine.

2] Versatile Color Selection:

Choose from 7 different colors

- | | |
|--|---|
| 1. Red  | 2. Blue  |
| 3. Yellow  | 4. Green  |
| 5. Red & Blue  | 6. Red & Yellow  |
| 7. Red, Blue, Yellow & Green  | |

Each color offers unique benefits.

3] Adjustable Intensity:

Customize the intensity level on a scale of 1 to 10, catering to your specific needs and preferences.

4] Pulse Selection:

Benefit from 11 pulse levels ranging from 10 to 110. Pulse functionality is available from 10 to 90, while 100 to 110 provides a steady, non-pulsating light.

5] Flexible Time Settings:

Set the treatment duration from 5 to 30 minutes across 6 different time intervals. The device automatically sets to 20 minutes upon startup.




biolight 7000sc

FEATURES

Safety Precautions

Read and follow all of the instructions in this manual before using. Proper use of the biolight 7000sc will prevent harm to the user and damage to private property. Only use the biolight 7000sc for the purposes specified in this manual.



Output and Color

First, it has a strong output that can enhance skin. The output is set differently for each color:



RED

40mW

Wavelength: 623nm



YELLOW

35mW

Wavelength: 590nm



GREEN

40mW

Wavelength: 512nm



BLUE

60mW

Wavelength: 470nm

You can choose from 7 different Options:

1. Red ●
2. Blue ●
3. Yellow ●
4. Green ●
5. Red + Blue ●●
6. Red + Yellow ●●
7. All: Red, Blue, Yellow & Green ●●●●

Intensity Selections

You can adjust the intensity from level 1 to 10.

► **PULSE OPTIONS:** Pulse can be adjusted in 11 levels from 10 to 110. For levels 10 to 90, it's pulse functionality. For levels 100 to 110, it's a steady light, and pulse functionality is 0.

- **10~30:** Muscle depth
- **40~60:** Subcutaneous layer
- **70~90:** Dermis layer
- **100~110:** Epidermal layer

Time Selections

Time can be adjusted in 6 levels from 5 to 30 minutes.

When the device is turned on, it's automatically set to 20 minutes.

- **MUSCLE RELAXATION:** 30 minutes
- **SUBCUTANEOUS LAYER:** 25 minutes
- **DERMIS LAYER:** 20 minutes

The theoretical basis of LED light therapy can be described as follows:

LED light therapy is based on the principle of photoactivation, where the application of LED light induces specific cellular responses without causing tissue damage. The therapy utilizes low-output LEDs that emit specific wavelengths of light, typically in the visible or near-infrared spectrum.

The interaction between the emitted light and the cells occurs through a process known as photobiomodulation. When the LED light is absorbed by chromophores within the cells, it triggers a series of biological reactions. One important chromophore is cytochrome c oxidase, which is part of the mitochondrial respiratory chain.

The absorbed light energy stimulates cytochrome c oxidase, leading to increased cellular energy production in the form of adenosine triphosphate (ATP). This heightened ATP production provides cells with the necessary energy to perform various functions, such as cellular repair, growth, and regeneration.





Basic biolight 7000sc Usage Instructions: biolight 7000sc with SF10 and SF3

- 1 Drink a glass of water (alkaline water is recommended).
 - 2 Have the customer lie down on the Biomat with the Acoustic-Mat on top.
(Acoustic Mat Settings: Relax 1, Bio-Mat Settings: 158°F).
 - 3 Wipe the face with a warm towel.
(Place a damp towel in the microwave and heat it for 40 seconds).
 - 4 Apply SF10 serum to the face in sequence and gently pat to absorb.
(Using the Facial Lifter device is recommended).
 - 5 Set the biolight 7000sc to intensity 10, color 1 (red), pulse 110, and a duration of 10 minutes.
(Use for 10 minutes after applying SF10 serum).
 - 6 After the 10-minute session, attach the SF3 Facial Mask to the face and connect the microcurrent battery.
 - 7 Set the biolight 7000sc to intensity 9, colors 7 (red, blue, yellow, and green), pulse 90, and a duration of 20 minutes.
(Use for 20 minutes after attaching the SF3 mask).
- By following these steps for a total of 30 minutes, Skin Care session will be completed.



SETTINGS FOR Rapid dermal Management


- ▶ **COLOR** Red, Blue, Yellow, and Green

 (Option 7 / 4 Colors)

- ▶ **INTENSITY** 90 mW

- ▶ **PULSE** 90 Hz

- ▶ **TIME** 20 Minutes





SETTINGS FOR Managing Outer Skin Layer

▶ COLOR	Red and Blue ●● (Option 5 / 2 Colors)
▶ INTENSITY	50 mW
▶ PULSE	100 Hz
▶ TIME	20 Minutes

1. RED
2. BLUE
3. YELLOW
4. GREEN
5. RED+BLUE
6. RED+YELLOW
7. ALL

Intensity of Light: 50 mW
Frequency Level: 100 Hz

TIMER

S Intensity Frequency ON/OFF

RICHWAY BIO LIGHT



SETTINGS FOR Whitening

- ▶ **COLOR** Blue
 ●
 (Option 2 / 1 Color)

- ▶ **INTENSITY** 80 mW

- ▶ **PULSE** 100 Hz

- ▶ **TIME** 20 Minutes

○ 1. RED

2. BLUE

○ 3. YELLOW

○ 4. GREEN

5. RED+BLUE

6. RED+YELLOW

7. ALL

TIME SET

20

▲

TIMER

Intensity
of Light

80

mW

Frequency
Level

100

Hz

S

▲

▲

⏻

Select
Color

Intensity

Frequency

ON
OFF

RICHWAY BIO LIGHT



SETTINGS FOR For Find Lines

- ▶ COLOR Red and Yellow
●●
(Option 6 / 2 Colors)

- ▶ INTENSITY 50 mW

- ▶ PULSE 110 Hz

- ▶ TIME 20 Minutes

● 1. RED

○ 2. BLUE

● 3. YELLOW

○ 4. GREEN

5. RED+BLUE

6. RED+YELLOW

7. ALL

S
Select Color

TIME SET

20

▲
TIMER

Intensity of Light

50 mW

Frequency Level

110 Hz

▲ Intensity ▲ Frequency

ON OFF

RICHWAY BIO LIGHT




SETTINGS FOR Resolving Skin Troubles

► COLOR Green
●
(Option 4 / 1 Color)

► INTENSITY 10 mW

► PULSE 110 Hz

► TIME 20 Minutes





SETTINGS FOR Treating Muscle Pain and Neuralgia

- ▶ COLOR Red
●
(Option 1 / 1 Color)

- ▶ INTENSITY 10 mW

- ▶ PULSE 40 Hz

- ▶ TIME 30 Minutes

1. RED

2. BLUE

3. YELLOW

4. GREEN

5. RED+BLUE

6. RED+YELLOW

7. ALL

TIME SET

30

▲
TIMER

Intensity
of Light

10

mW

Frequency
Level

40

Hz

S

▲

▲

⏻

ON
OFF

RICHWAY BIO LIGHT

PREMIUM FACE SERUM

SF10™ *with* Penellagen

BIOTECHNOLOGY

A Quality Life is not only about living longer,
but also living a healthier life.

With Biotechnology,
we can open the doors of challenge
the reversing of aging.

53 Ingredients



MICRO CURRENT

SF3™ Facial Mask

Microcurrent stimulates the face by
sending gentle electrical waves
through the skin and tissues to
the facial muscles.

Microcurrent has been shown to
stimulate ATP production,
stimulate facial muscle vitality
and promote collagen and
elastin production.

The microcurrent pair stimulates
muscle movement and accelerates the
synthesis of collagen on the face.

